




### The Universal Solution for PICC Securement

The SorbaView SHIELD Plus system is an effective, easy to use method for securing all PICC catheters – *including silicone*. SHIELD Plus is a great solution for hospitals looking to standardize PICC dressing change protocols and care. SorbaView SHIELD Plus is one more example of Centurion's customer-centered approach to product development.


## BENEFITS FOR CAREGIVERS, HOSPITALS AND PATIENTS

- Effective, easy-to-use securement system for all PICCs – prevents catheter movement
- Works with **silicone** catheters and polyurethane catheters
- System is easy to use as well as economical
- 7-day wear time means fewer dressing changes on central lines
- Absorbent pad wicks away moisture
- Combines the benefits of SorbaView SHIELD's lasting dwell times and patient-friendly adhesive with HubGuard's enhanced catheter securement
- Provides cushion, increased breathability and patient comfort
- Establishes a standard of PICC care hospital wide
- Contoured design of the HubGuard works well with antimicrobial disks
- System designed to help lower the risk of catheter-related bloodstream infections (CRBSIs) associated with central lines

## HUBGUARD APPLICATION INSTRUCTIONS



## SORBAVIEW SHIELD APPLICATION INSTRUCTIONS



**CENTURION**  
MEDICAL PRODUCTS

100 CENTURION WAY | P.O. BOX 510 | WILLIAMSTON | MICHIGAN | 48895  
517.546.5400 | 800.248.4058 | fax 517.546.9388 | www.centurionmp.com

©2010 Centurion Medical Products Corporation LIT131V1

Product #SV353UDTHG  
• 25 per subpack  
• 50 per case